



QPRT Rabbit pAb

Catalog No	YP-Ab-19466
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	QPRT
Immunogen	Synthesized peptide derived from human QPRT. AA range:119-219
Specificity	This antibody detects endogenous levels of QPRT at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	QPRT;Nicotinate-nucleotide pyrophosphorylase [carboxylating] (Quinolate phosphoribosyltransferase [decarboxylating]) (QAPRTase) (QPRTase)
Observed Band	
Calculated Molecular Weight	33kD
Cell Pathway	
Tissue Specificity	
Function	Involved in the catabolism of quinolinic acid (QA) .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

