



CAP2 Rabbit pAb

Catalog No	YP-Ab-19550
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	CAP2
Immunogen	Synthesized peptide derived from human CAP2. AA range:22-122
Specificity	This antibody detects endogenous levels of CAP2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	CAP2;Adenylyl cyclase-associated protein 2 (CAP 2)
Observed Band	
Calculated Molecular Weight	52kD
Cell Pathway	Cell membrane ; Peripheral membrane protein .
Tissue Specificity	
Function	May have a regulatory bifunctional role.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

