



TTYH1 Rabbit pAb

Catalog No	YP-Ab-19551
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TTYH1
Immunogen	Synthesized peptide derived from human TTYH1. AA range:50-120
Specificity	This antibody detects endogenous levels of TTYH1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	TTYH1;Protein tweety homolog 1 (hTTY1)
Observed Band	
Calculated Molecular Weight	50kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in brain , eye , ovary and testis , and at lower levels in muscle , placenta , liver and lung.
Function	Probable chloride channel. May be involved in cell adhesion (By similarity) . ; Isoform 3 may be a Ca (2+) -independent and swelling-activated chloride channel , possibly involved in regulation of cell volume.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images