



CALHM2 Rabbit pAb

Catalog No	YP-Ab-19558
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	CALHM2 FAM26B
Immunogen	Synthesized peptide derived from human CALHM2. AA range:182-252
Specificity	This antibody detects endogenous levels of CALHM2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	CALHM2 FAM26B;Calcium homeostasis modulator protein 2 (Protein FAM26B)
Observed Band	
Calculated Molecular Weight	36kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	Pore-forming subunit of a voltage-gated ion channel.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

