



FUT4 Mouse mAb

Catalog No	YP-mAb-19545
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive .
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	LeX; CD15; ELFT; FCT3A; FUTIV; SSEA-1; FUC-TIV; FUT4
Observed Band	59kDa
Calculated Molecular Weight	59kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	The product of this gene transfers fucose to N-acetyllactosamine polysaccharides to generate fucosylated carbohydrate structures . It catalyzes the synthesis of the non sialylated antigen, Lewis x (CD15).
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

