



FAH Mouse mAb

Catalog No	YP-mAb-19552
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	FAH
Immunogen	Synthesized peptide derived from human FAH. AA range:229-329
Specificity	This antibody detects endogenous levels of FAH at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	FAH;Fumarylacetoacetase (FAA) (Beta-diketonase) (Fumarylacetoacetate hydrolase)
Observed Band	
Calculated Molecular Weight	46kD
Cell Pathway	
Tissue Specificity	Mainly expressed in liver and kidney. Lower levels are also detected in many other tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

