




NOP2 Mouse mAb

Catalog No	YP-mAb-19553
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive .
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	NOL1; p120; NSUN1; NOP120; NOP2
Observed Band	110kDa
Calculated Molecular Weight	89kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables RNA binding activity . Involved in positive regulation of cell population proliferation; regulation of signal transduction by p53 class mediator; and ribosomal large subunit assembly . Located in nucleolus .
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of various lysates using NOP2 Mouse mAb at 1:1000 dilution .

Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) at 1:10000 dilution .

Lysates/proteins: 25μg per lane .

Blocking buffer: 3% nonfat dry milk in TBST .

Detection: ECL Basic Kit (RM00020).

Exposure time: 30s .