



PADI2 Rabbit pAb

Catalog No	YP-Ab-18735
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	PAD2; PDI2; PADI2
Immunogen	Synthesized peptide derived from human PADI2
Specificity	This antibody detects endogenous levels of PADI2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PAD2;PDI2;PAD-H19;PADI2;Protein-arginine deiminase type-2 (PAD-H19) (Peptidylarginine deiminase II) (Protein-arginine deiminase type II)
Observed Band	75kDa
Calculated Molecular Weight	76kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Detected in keratinocytes in epidermis (at protein level).
Function	Catalyzes the deimination of arginine residues of proteins.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

