



RFA1 mouse mAb

Catalog No	YP-mAb-12518
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPA1 REPA1 RPA70
Immunogen	Synthesized peptide derived from human RFA1
Specificity	This antibody detects endogenous levels of Human,Monkey,Mouse,Rat RFA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSSB;RF-A;RP-A;REPA1;RPA70;MST075;PFBMFT6;RPA70/RPA1;RFA1; RPA1 REPA1 RPA70;Replication protein A 70 kDa DNA-binding subunit (RP-A p70;Replication factor A protein 1;RF-A protein 1;Single-stranded DNA-binding protein)
Observed Band	
Calculated Molecular Weight	68kD
Cell Pathway	Nucleus . Nucleus, PML body . Enriched in PML bodies in cells displaying alternative lengthening of their telomeres. .
Tissue Specificity	
Function	function:It participates in a very early step in initiation. RP-A is a single-stranded DNA-binding protein. Absolutely required for simian virus 40 DNA replication in vitro.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the replication factor A protein 1 family., subunit:Heterotrimer of 70, 32 and 14 kDa chains. The DNA-binding activity may reside exclusively on the 70 kDa subunit. Interacts with RIP and XPA. Interacts with RPA4.,



Background

function:It participates in a very early step in initiation. RP-A is a single-stranded DNA-binding protein. Absolutely required for simian virus 40 DNA replication in vitro.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the replication factor A protein 1 family., subunit:Heterotrimer of 70, 32 and 14 kDa chains. The DNA-binding activity may reside exclusively on the 70 kDa subunit. Interacts with RIP and XPA. Interacts with RPA4.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

